

Sample courier guide

Please use the thermocol or cardboard box with ice pack. Fill the test request form and choose the tests you require for your patient.

To ensure good quality samples

Collect the samples as directed in test request form. Test tubes and amount of sample are indicated on the test request form against the required test. Next label the tube with patient name, Age, gender.

Samples should be collected late evening or prior to dialysis (if patient on dialysis) and given to courier service to be flown out early morning. This will ensure that the samples reach us by the next day afternoon in good condition.

To be couriered to the following address

ATDI Lab Pvt Ltd
6 B, 6th Floor, Trust House 35,
Dr. E. Borges Road Opp Shirodkar High School,
Parel, Mumbai 400012

The form and related documents can be sent along with sample and a scanned copy soft copy on email.

NEFT Details for online payment

Bank account No- 6645209310
IFSC code- KKBK0001388
MICR Code- 400485106
Branch- Chembur, Mumbai

Contact us on

Email atdi5353@gmail.com, sjogaleatdi@gmail.com
Lab Manager, Ms Laxmi Chowti – 9321757662